

PHLDB1-pleckstrin homology-like domain, family B, member 1 Gene

Reference: PTXBC014510

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pleckstrin homology-like domain family B member 1; LL5alpha; pleckstrin homology-like domain family B member 1 variant 3; pleckstrin homology-like domain family B member 1 variant 4; protein LL5-alpha; pleckstrin homology like domain family B member 1		

For research use only.

GENE DESCRIPTION	pleckstrin homology-like domain, family B, member 1
GENE ID	23187
ORIGIN SPECIES	Human
ACCESSIONS	BC014510
VECTOR	pDONR223
NCBI SYMBOL	PHLDB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgcaagattaaatcatggaagaagcgtgtttgtcttcgaccggctcaagcgcacccttcctattatgtgacaagcatgagacgaagctgaaggagtcattctattccaggccattgaggaagtgtactacgaccacctgcgcagtgcaagcaagagccgaaccagccc
tcaccttcgctaaagaccatgaccggctgtactacatggtggcccatctgcagaggccatgcgtatctggatggatgcatgtcagagggtgagggtactacactcagttcatgaactaa