

MB-myoglobin Gene

Reference: PTXBC014547

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PVALB; myoglobin		
For research use only.			

GENE DESCRIPTION	myoglobin
GENE ID	4151
ORIGIN SPECIES	Human
ACCESSIONS	BC014547
VECTOR	pDONR223
NCBI SYMBOL	MB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctcagcgacggggaatggcagttggtgctgaacgtctgggggaaggfagggtgacatccaggccatgggcaggaaagtcctcaggtctttaaagggtcaccagagactctggagaagtttgacaagttcaagcacctgaagtcagaggacgagatgaaggcatctg
aggacttaagaagcatggtgccactgtgctcaccgcccgggtggcatccttaagaagaagggtcatcatgaggcagagattaaagccctggcacagtcgcatgccaaagcacaagatccccgtgaagtacctggagttcattcgaatgcatcatccaggttctgcagagcaag
catccccgggactttggtgctgatccaggggccatgaacaaggccctggagctgttccggaaggacatggcctccaactacaaggagctgggcttccagggttag