



Reference: PTXBC014566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDC15; disintegrin and metalloproteinase domain-containing protein 15; MDC-15; a disintegrin and metalloproteinase domain 15 (metargidin); metalloprotease RGD disintegrin protein; metalloproteinase-like, disintegrin-like, and cysteine-rich protein 15; ADA		

For research use only.

GENE DESCRIPTION	ADAM metallopeptidase domain 15
GENE ID	8751
ORIGIN SPECIES	Human
ACCESSIONS	BC014566
VECTOR	pDONR223
NCBI SYMBOL	ADAM15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]