

## C15orf15-chromosome 15 open reading frame 15 Gene

Reference: PTXBC014576

|                |                                                                                                                                                                                                                                    |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C15orf15; HRP-L30-iso; L30; RLP24; RPL24; RPL24L; TVAS3; 60S ribosomal protein L30 isolog; homolog of yeast ribosomal like protein 24; my024 protein; ribosomal L24 domain-containing protein 1; ribosomal L24 domain containing 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 15 open reading frame 15                                                               |
| <b>GENE ID</b>          | 51187                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC014576                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C15orf15                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgcgtatcgaaaagtgttattctgttcggggcccatctatctggacacggcatgatgttcgtccgcaacgattgcaagggttcggtttgcaaatctaaatgtcataaaactttaaaaagaagcgcaatcctcgcaagtttaggtggaccaaaagcattccggaaagcagctgtaaaagagcttacagtggataatcatttgaatttgaaaaacgtagaatgagcctatcaataaccagcgagagctatggaataaaactattgatgcgatgaagagagttgaagaaatcaaacagaagcgccaagctaaattataatgaacagattgaagaaaaataaagagctacagaaagttcaggat atcaaaagaagtcagcaaaacatccatcttatccgagccctcttgcaggcaagggaacagttggaagagaaaaatggtacagcagttacaagaggatgtggacatggaagatgctccttaa