

CEP170L-centrosomal protein 170kDa-like Gene

Reference: PTXBC014590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	centrosomal protein 170kDa-like
GENE ID	645455
ORIGIN SPECIES	Human
ACCESSIONS	BC014590
VECTOR	pDONR223
NCBI SYMBOL	CEP170L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctactagtcttcattcaaacaccggattaaagcagggaagactacatccgagattggactgctcatcgagaagagatagccaggatcagccaagatcttgcctcattgctcgggagatcaacgatgtagcaggagagatagattcagtgacttcacaggcactgccctagtacc
acattggtgatcgtgtttttatgaagcctcaacttccaaagattcctcattgctcattccaaacaccagaagggaacaacggctgctggtgatccaagacctcaagcagcagagcctccgatcacttaacaattacaaggcggagaacctggagcaggggatgaagtcattggga
gataatctgctgctgcatccgtctttcagttccttaagaagataagacaatctatagataagacagctggaagatcagaatatttttaagacaaagatcggaattgggatgacatagaagcaaatgaagagccgaaagtgaagtcctattgtgaaacctcagcatggagatttctta
tcttacaggaaactgaaagagtagaaaagcagctacaagcaatcaatgctatgattgatcctgatggaacttggaggctctgaacaacatggatttccagtgctatgttgccatcaccgaaacagaagtcagccctgtgataaacaccacagccgggtcagacaccaacacttg
gccaaccagaagctagggtcttcattcctgctgctgtttacggcagctgaattgagaatgctgaatctgaggtgatttcagtatacttcaatagatcaacctgatggggaagaggaagatgtacagtacataaatga