



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SMCP-sperm mitochondria-associated cysteine-rich protein Gene

Reference: PTXBC014593

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSMCSGEN1; MCS; MCSP; sperm mitochondrial-associated cysteine-rich protein; mitochondrial capsule selenoprotein; testicular tissue protein Li 119; testicular tissue protein Li 161; sperm mitochondria associated cysteine rich protein		

For research use only.

GENE DESCRIPTION	sperm mitochondria-associated cysteine-rich protein
GENE ID	4184
ORIGIN SPECIES	Human
ACCESSIONS	BC014593
VECTOR	pDONR223
NCBI SYMBOL	SMCP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgaccagacaaaacacagtaaatgctgcccagcaaaaggcaatcaatgctgccaccacagcagaaccagtgctgcccagcaaaaaggcaatcaatgctgccaccaaaacagaaccagtgctgcccagcaaaaaggcagtaatgctgccaccaaaacacaaatcactgctgcc
agccaaaacccccafgctgctcattcagccgagtgctgtgtgtttgagaccaagcctgaagtctcaccccttaacatgagtgctgagcccaactcaccgcaactcaggacaaagggctgtaaccaaccagcagcagcccccatagcccaaaatgagtccagcccaagcaaatga