

C7orf34-chromosome 7 open reading frame 34 Gene

Reference: PTXBC014596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf34; ctm-1; MSSP-binding protein CTM-1; chromosome 7 open reading frame 34		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 34
GENE ID	135927
ORIGIN SPECIES	Human
ACCESSIONS	BC014596
VECTOR	pDONR223
NCBI SYMBOL	C7orf34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctccctggccccagctctgcagggcagtggtctctgttctctatcttgcctgctgcaggtgaagcctctgaacgggagccagggcccaagatgggagccagacagagaaaacgcctctgcagaccagaatcaagaacagttcgaagagcactttgtggcctcctcagtg
ggtagatgtggcaggtgggtgacatggccagcaggaagaagaccagtcgtccaagacggcagctgttcacaagcactcttccacctcagctctgttagtctggccagtgctcatggtttctcaggaggccatfgaggcggaatcccaaatccaactctgttcatgctcactc
actga