

GLIPR1L1-GLI pathogenesis-related 1 like 1 Gene

Reference: PTXBC014603

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALKN2972; PRO7434; GLIPR1-like protein 1; GLI pathogenesis related 1 like 1		
For research use only.			

GENE DESCRIPTION	GLI pathogenesis-related 1 like 1
GENE ID	256710
ORIGIN SPECIES	Human
ACCESSIONS	BC014603
VECTOR	pDONR223
NCBI SYMBOL	GLIPR1L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctgaagaataattcagttgtttatggatcttgggtctgtgtttgtagccactacatctcacaatcccatccatcactgacccacactttatagacaactgcatagaagcccaacgaatggcgtggcaagtcacccctcccgccgacatgaatacatgattgggataaa
ggtttagcaaatggtctaaagcatgggcaaacagtgcaatttgaacataatgactgtttggataaatcatataaatgctatgcagcttttgaatatgttggagaaaatctggttaggtggaataaagtcattcacaccaagacatgccattacggcttggataatgaacccaattttatgat
tttgatagctctatcatgctccagagctgtgtggccattatcacagtttagttgggccaattcattttatgtcgggtgtgcagttgcaatgtgtcctaacttggggagcttcaactgcaatatttggatgcaactacggacctgcaggaaatttgcataatgcctcttacgtaagaggagaatcttg
ctctctctgctcaaaagaagagaaatgtgtaagaacctctgcataaatccatttctgaagccaacggggagagcacctcagcagacagccttaaccattcagcttaggttttcttctctgagaatcttttaa