



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C19orf36-chromosome 19 open reading frame 36 Gene

Reference: PTXBC014609

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf36; IMAGE:4215339; izumo sperm-egg fusion protein 4; sperm 22 kDa protein c113; IZUMO family member 4		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 36
GENE ID	113177
ORIGIN SPECIES	Human
ACCESSIONS	BC014609
VECTOR	pDONR223
NCBI SYMBOL	C19orf36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]