

PDPN-podoplanin Gene

Reference: PTXBC014668

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGGRUS; GP36; GP40; Gp38; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; T11A; PA2.26 antigen; T1-alpha; glycoprotein 36; glycoprotein, 36-KD; hT1alpha-1; hT1alpha-2; lung type I cell membrane associated glycoprotein; lung type-I cell membrane-associated glycopr		

For research use only.

GENE DESCRIPTION	podoplanin
GENE ID	10630
ORIGIN SPECIES	Human
ACCESSIONS	BC014668
VECTOR	pDONR223
NCBI SYMBOL	PDPN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaaggtgtcagctctgctctctgtttgggaagcgcgtcgtctgggtctggcagaaggagccagcacagccagccagaagatgacactgagactacaggtttggaagcggtgtgcatgccaggtgccgaagatgatgtgtgactccagggaaccagcgaagaccgct
ataagctctggttgacaactctggtggcaacaagtgtcaacaggtgaacaggcattcgcatcgagatctgccaacttcagaagcacagtcacgcgcaagaacaagtcgaagcgcacagcctcaaacgtggccaccagtcactccacggagaagtggtggagacacagac
aacagttgagaagatggtttgtcaacagtgacctggttggaatcatagttgggtcttactagccatcggtctcattgtggaatcatcgtttgtgttatcgaaaaatgtcgggaaggtactgcctaa