

MANBAL-mannosidase, beta A, lysosomal-like Gene

Reference: PTXBC014672

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein MANBAL; mannosidase beta A like; mannosidase, beta A, lysosomal-like; mannosidase beta like		

For research use only.

GENE DESCRIPTION	mannosidase, beta A, lysosomal-like
GENE ID	63905
ORIGIN SPECIES	Human
ACCESSIONS	BC014672
VECTOR	pDONR223
NCBI SYMBOL	MANBAL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctctgacctagacttctcacctccggaggtgcccagagcccactttctggagaacctgtacgggtacgactcttctgggagccattctccagctcatctgtgtgtgtggtccatcatgtaccattcccaagtcacagagcgaggctgaacctgtgagcccagaagtgtga
ggtagcagaggaagcccaaggctgtcttctgtgaacaagagcccaagaaagagactaagaagaagcggtag