



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PHB2-prohibitin 2 Gene

Reference: PTXBC014766

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

BAP; BCAP37; Bap37; PNAS-141; REA; hBAP; p22; prohibitin-2; B-cell associated protein; B-cell receptor-associated protein BAP37; D-prohibitin; repressor of estrogen receptor activity; prohibitin 2

For research use only.

GENE DESCRIPTION	prohibitin 2
GENE ID	11331
ORIGIN SPECIES	Human
ACCESSIONS	BC014766
VECTOR	pDONR223
NCBI SYMBOL	PHB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgagaactggaaggacttggcggggcggctgccgcggggcccccgggcacgtggcgccgctgaagctgttgcctggggcggcgccgtggcctacgggtgtgcggaactgtgttaccgttggaaaggcgggcacagagccatctcttcaatggatcgtgtgagtgca
caggacacatctctggccggagcctctacatcagatcgtgttgcacacccatctatctgaatcttggcggcagacacgtgaaatactctctccctacagacgtcgaagacacacagatgttgaatcttccctcgaggtgtgttctccgacccaatctcaggaatctctcagatct
accagccgttagggctgactacaggaagcaggtgttgcctgtcttcaacgaggtgtccaaggtgtgtgttggggccaaagtcaatccctacacagcaccgacggcgccgacgtatctctgtctacggcggaatctcagagagagggcccaagatcagctcatctcgaat
gatgtggccatcacagagctgagctttagccgagagtacacagctgctgttagaagcccaacaagtggcccagcaggaggccggcggcccaattctgttgaataaaagcaagcagggaacagcggcgaaaattgtcaggccggagggtgagccggaggctgccaagatctgt
gagaagcactgagcaagaacctggctacatcaaatctcgaaagtctgagcagccagaatatctccaagcagatccacatcacagaatctgtatctatctacagctgacaacctgtgtcgtgaacctacagatgaagatttcaccaggggaagtgacagctctcaaggttaagaa
ata

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1