

SLC25A6-solute carrier family 25 (mitochondrial carrier; adenine nucleotide translocator), member 6 Gene

Reference: PTXBC014775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AAC3; ANT; ANT 2; ANT 3; ANT3; ANT3Y; ADP/ATP translocase 3; ADP,ATP carrier protein; ADP,ATP carrier protein 3; ADP,ATP carrier protein, liver; ADP/ATP translocator of liver; adenine nucleotide translocator 3; solute carrier family 25 (mitochondrial carr		

For research use only.

GENE DESCRIPTION	solute carrier family 25 (mitochondrial carrier; adenine nucleotide translocator), member 6
GENE ID	293
ORIGIN SPECIES	Human
ACCESSIONS	BC014775
VECTOR	pDONR223
NCBI SYMBOL	SLC25A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggaacagggccatctcctcgccaaagactcttgccggaggcgcgcgcgcacatctccaagacggcgtggctccgacgagcgggtcaagctgctgctgcaagtcacgacgcccagcagatcgccgcccagaagcagtagaaggcgcgtgactgcatgtgc
cgcatcccaaggagcagggcgctgctctctgaggggcaaccttgccaacgtcattcgtacttcccactcaagccctcaacttcgcttcaaggataagtagaagcagatctctctggggcgctggacaagcacacgagttctggaggtaattgcccggcaacctggcctccg
gggtgcggcgccgacgtccctctgcttctgtacccgctggattttccagaacccgctggcagcggagctgggaagtcaggcacagagcgcgagttccgagggcctgggagactgcctggtgaagatcaccaagtcgacggcatccgggctgtaccagggttcagtg
tctccgtgagggcgcacatcatctaccggcgccctacttggcgtgtacgatacggccaaggcgcgtcccccagaccccaagaacacgcacatcgtggtgagctggtgatgcgcgacacgtgacggcgtggccggtggtgtctacccttcacacgggtgcggcgccgcat
gatgatgcagttcgggcgcaaggagctgacatcatgtacacgggcaccgtcgtgttgagggaagatcttcagagatgagggggcgaaggccttctcaagggtgctgtgtccacgtcctcggggcagtgggggcgcccttcgtgtgtcctgtacgacgagctcaagaagggtg
atctaa