

LOC541471-hypothetical LOC541471 Gene

Reference: PTXBC014776

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGD2; LINC00978; MIR4435-1 host gene (non-protein coding)		
For research use only.			

GENE DESCRIPTION	hypothetical LOC541471
GENE ID	541471
ORIGIN SPECIES	Human
ACCESSIONS	BC014776
VECTOR	pDONR223
NCBI SYMBOL	LOC541471
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagacactgaaatcacgactcatccccctcagcacctctactgttgcgccgatcacagccggaatgcagctgaaagattccctggggcctggttccaactgccactgtggactctgaggcctctgcatttgcgggtggtctgcctgtga