

C9orf80-chromosome 9 open reading frame 80 Gene

Reference: PTXBC014881

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf80; HSPC043; MISE; SOSSC; SSBIP1; hSSBIP1; SOSS complex subunit C; SSB-interacting protein 1; hSSB-interacting protein 1; minute INTS3/hSSB-associated element; sensor of single-strand DNA complex subunit C; sensor of ssDNA subunit C; single-stranded		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 80
GENE ID	58493
ORIGIN SPECIES	Human
ACCESSIONS	BC014881
VECTOR	pDONR223
NCBI SYMBOL	C9orf80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcaactcttcaggacaaggtttcaaaacaaaaatagagttgcaatcttggcagaactggacaagagagaaaaaaaactactatgcagaaccagtcttcaacaatcatcctggagctagcattgcactctcgagaccctctcttaataaggacttcgggatcacgtgagcag
cagcatattgcagcccaacagaaggcagctttgcagcatgctcatgcacattcatctggatacttcatcactcaagactctgcatctgggaaccttattcttctgttttacctcgcttgaccagaatga