



Reference: PTXBC014895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	aldehyde dehydrogenase family 16 member A1; aldehyde dehydrogenase 16 family member A1		

For research use only.

GENE DESCRIPTION	aldehyde dehydrogenase 16 family, member A1
GENE ID	126133
ORIGIN SPECIES	Human
ACCESSIONS	BC014895
VECTOR	pDONR223
NCBI SYMBOL	ALDH16A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]