



Reference: PTXBC014904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LRDD; PIDD; p53-induced death domain-containing protein 1; leucine-rich repeat and death domain-containing protein; leucine-rich repeats and death domain containing; p53-induced protein with a death domain; p53-induced death domain protein 1		

For research use only.

GENE DESCRIPTION	leucine-rich repeats and death domain containing
GENE ID	55367
ORIGIN SPECIES	Human
ACCESSIONS	BC014904
VECTOR	pDONR223
NCBI SYMBOL	LRDD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]