



Reference: PTXBC014921

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SFRS15; SRA4; splicing factor, arginine/serine-rich 15; CTD-binding SR-like protein rA4; SR-like CTD-associated factor 4; pre-mRNA splicing SR protein rA4; splicing factor serine alanine 15; SR-related CTD associated factor 4		

For research use only.

GENE DESCRIPTION	splicing factor, arginine/serine-rich 15
GENE ID	57466
ORIGIN SPECIES	Human
ACCESSIONS	BC014921
VECTOR	pDONR223
NCBI SYMBOL	SFRS15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr