



Reference: PTXBC014945

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIMAH2; armadillo repeat-containing protein 5; armadillo repeat containing 5		
For research use only.			

GENE DESCRIPTION	armadillo repeat containing 5
GENE ID	79798
ORIGIN SPECIES	Human
ACCESSIONS	BC014945
VECTOR	pDONR223
NCBI SYMBOL	ARMC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]