

C20orf107-chromosome 20 open reading frame 107 Gene

Reference: PTXBC014951

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf107; C20orf107; dJ1153D9.4; protein FAM209B; family with sequence similarity 209 member B		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 107
GENE ID	388799
ORIGIN SPECIES	Human
ACCESSIONS	BC014951
VECTOR	pDONR223
NCBI SYMBOL	C20orf107
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggacgctgaaatcgtccctggtcctgcttctgtgctcactgcagctatgctttatgttctctctctgagacagaaaactagcgaacccaggagggtgccgtgtggagagcacttctggattcggcagaacctaccagagcacaccaaggctggcttgggagcaaatggctc
tggcttttggctgtgtgctccgtttgtgatactgaagtgtcaagagacagtgagaagaataaggaaggatggctcatttttttacaccatattgattcaatctcaggagctcagggaacggatgtcttag