

## C6orf115-chromosome 6 open reading frame 115 Gene

Reference: PTXBC014953

|                |                                                                                                                        |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C6orf115; Costars; HSPC280; PRO2013; costars family protein ABRACL; ABRA C-terminal-like protein; ABRA C-terminal like |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 6 open reading frame 115                                                               |
| <b>GENE ID</b>          | 58527                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC014953                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C6orf115                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaatgtggatcagaggttaacctcttagtggaggaaattcatgtttgggtcaaaaaatgctgatggaaagtaagcgtgaaattggggctctctccgtgatgataaatgtgccaacctcttgaagcattgtaggaactcttaagctgcaaacgaaggagattgaacatccag  
gagagctgctctgcaaggtgttcacatgatgttgacattatattactgcaagattaa