

## SNAPC1-small nuclear RNA activating complex, polypeptide 1, 43kDa Gene

Reference: PTXBC014984

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PTFgamma; SNAP43; snRNA-activating protein complex subunit 1; PSE-binding factor subunit gamma; PTF subunit gamma; SNAPc 43 kDa subunit; SNAPc subunit 1; proximal sequence element-binding transcription factor subunit gamma; small nuclear RNA activating co |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | small nuclear RNA activating complex, polypeptide 1, 43kDa                                       |
| <b>GENE ID</b>          | 6617                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC014984                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SNAPC1                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggggactctcccgccctgcagaccgactgcgagggcgtgctcagccgcttcaggagacggagcagtgctacgttcgaggacttcacggagctctggagaacatgaagttcgggactatcttctgtggcagaatgagaaatttagaaaaaacatgtttcaaaaagaagcttagcttt  
ggcttggcgatatattttacctccatacacttccagatcagagttgggtgctttgctatctgctatggtattatataataacccaactgtgtcaacccaaacaaaagatcagagttgccctgaagattgggatgaagttttaaaatttcagcaagatttagtaaatgcacagcatttgcagcttataat  
tttaggaagctacgactagacagagcatttcacitacagcaatgccaaattgctgctataggtgagaaaaaaattcaccgagctgaagttacagaagaatttaaggaccacagtgatcgtgtgatgaaacttaccctctgattatagaggaatgctgaatgtcatgatcattatc  
agaacatgaacatgtaatttcagttgataagtcgaagccagataaagccctcagcttgataaaggtgattttttgacaattatgaacatagtttggagcatcagcagtgccacaaagacagaagaatccatccttaagtcataaactaatgatggagaagaaaaatggaaggaaat  
tcacaagaacggagagatgtgaaggcgaatcattagcgaataaaatcaaggccttttcagttgtacacagcgcataatcaagaaggcctcgaagtcacactcgactctctgactctgattctgcatctgggtcaaggcgaagtcgaagcaactaggaagaaagagaagaa  
agaagattgaaccagcaggaaggaagatgtctcagaaacaaggcaatgtgcagaatatacacaaggagataaaccttaagtctgagtgctgtattacagaagaagaagaatgaaagtttgagtggaacagagttcactgcacccaagaaggagaaaacactga