

RAB18-RAB18, member RAS oncogene family Gene

Reference: PTXBC015014

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB18, member RAS oncogene family; RAB18 small GTPase; RAB18LI1; WARBM3; ras-related protein Rab-18		

For research use only.

GENE DESCRIPTION	RAB18, member RAS oncogene family
GENE ID	22931
ORIGIN SPECIES	Human
ACCESSIONS	BC015014
VECTOR	pDONR223
NCBI SYMBOL	RAB18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgaggacgtgctaaccacccctgaagatcctcatcgcgcgagagtggtgggcaagtcagcctgctcttgaggttcacagatgatacgtttgatccagaactgcagcaacaataggtgtgactttaaggtgaaacaatttcagtgatgaaataaggctaaactgcaatat
gggatactgctggtcaagagaggtttagaacatttaactccagctattatagaggtgcacaggggtgtatattagttatgatgtcacaagaagagatactttgtaaaactggataattgggttaaatgaattggaacatactgtacaagaatgacatagtaaacatgctagtggaaataaaatc
gataaggaaaatcgtgaagtcgatagaatgaaggcctgaaatttgcacgaaagcattcattgtatttatagaggcaagtgcaaaaactgtgatggtgtacaatgtgcctttgaagaactgttgaaaagatcattcagaccctggactgtgggaaagtgaaccagaataaaggagtc
aaactgtcacacaggggaagccaaggaggaggagcctgtggtgtatttgcctctgtgtataa