

ATP5L-ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit G Gene

Reference: PTXBC015128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP5JG; ATP synthase subunit g, mitochondrial; ATP synthase g chain, mitochondrial; ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit G; ATP synthase, H ⁺ transporting, mitochondrial F1F0, subunit g; ATPase subunit G; F1F0-type ATP synthase		

For research use only.

GENE DESCRIPTION	ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit G
GENE ID	10632
ORIGIN SPECIES	Human
ACCESSIONS	BC015128
VECTOR	pDONR223
NCBI SYMBOL	ATP5L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaattgtccgtaacctgtggagaagaccccgccgtggtgaacgtgctgtgacttactcgaagcctcgattggccacatttggctactacgccaaggttgagctggttcctccaccctgctgagatccctagagctattcagagcctgaaaaaataagccaatagtgtctcaga
ctggtagcttcaaacagctcacagttaaggaagctgtgctgaatggtttggtggccactgaggtgttgatgtggtttatgtcggagagattataggcaagcggggcatcattggctatgatgttga