

CYCS-cytochrome c, somatic Gene

Reference: PTXBC015130

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYC; HCS; THC4; cytochrome c; cytochrome c, somatic		

For research use only.

GENE DESCRIPTION	cytochrome c, somatic
GENE ID	54205
ORIGIN SPECIES	Human
ACCESSIONS	BC015130
VECTOR	pDONR223
NCBI SYMBOL	CYCS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgatgttgagaaaggcaagaattttattatgaagtgttccagtaaccacacgttgaaggaggagcaagcacaagactgggccaatctccatgggtctctttggcggaagacaggtcaggccctggatactttacagccgccaataagaacaaggcatcatctggg
gagaggatacactgatggagtatttgagaatccaagaagtacatccctggaacaaaaatgatcttgcggcattagaagaaggagaaagggcagacttaatagtctatctcaaaaagctactaatgagtaa