



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

EXOC7-exocyst complex component 7 Gene

Reference: PTXBC015165

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

2-5-3p; EX070; EXO70; EXOC1; Exo70p; YJL085W; exocyst complex component 7; exocyst complex component Exo70

For research use only.

GENE DESCRIPTION	exocyst complex component 7
GENE ID	23265
ORIGIN SPECIES	Human
ACCESSIONS	BC015165
VECTOR	pDONR223
NCBI SYMBOL	EXOC7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1