

ALKBH8-alkB, alkylation repair homolog 8 (E. coli) Gene

Reference: PTXBC015183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABH8; TRM9; TRMT9; alkylated DNA repair protein alkB homolog 8; AlkB homologue 8; S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8; alkB, alkylation repair homolog 8; tRNA (carboxymethyluridine(34)-5-O)-methyltransferase ABH8; tRNA methyltran		

For research use only.

GENE DESCRIPTION	alkB, alkylation repair homolog 8 (E. coli)
GENE ID	91801
ORIGIN SPECIES	Human
ACCESSIONS	BC015183
VECTOR	pDONR223
NCBI SYMBOL	ALKBH8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagcaaccatcaagtaattacaactcagtaaaactgagaagaagttcttaaggaaacagattaaagccaagcactcttctgagacatgaaggcattgagacagatcctatgccactcagagcctggtgtgtccaatggtgtgttggtggtatggtgtgagtcgggaaccagctgc
tccccgttttagagaaatgtggactggtggatgctctcttaattgccacctaacaagccgtactcatttgaagatacagaactacagaagaatctaagagagcctatgttacccctcaatggaaaaagtaagtggatgatttaggacaaaagatcactctgtatttgaatttggtaaaagtgc
gtggaaggagttgaggcctcaagcctaccaggactcatggtagtagaagaataattctctgaggaggagaaaatgcttttggaaagtgtgattggacagaagatacagacaatcaaaactctcaaaatccttaaaacacagagagtaagcattttggtatgagttccactatg
agaacaacaatgtagataaagataagccattatctgggggtctctgacatttggaaagccttttggagaatggttgaggaaaggttacattaaacataaacctgatcaaatgaccataatcagtatgaacctgggcaagattgcatgattttaa