

HHIPL2-HHIP-like 2 Gene

Reference: PTXBC015201

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1822L; HHIP-like protein 2; hedgehog interacting protein-like 2; HHIP like 2		
For research use only.			

GENE DESCRIPTION	HHIP-like 2
GENE ID	79802
ORIGIN SPECIES	Human
ACCESSIONS	BC015201
VECTOR	pDONR223
NCBI SYMBOL	HHIPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggcccaagatcctcagggtgctgcagctgcctgagcgaggtggccaacgggctgaggaaacccgtctccatggtccatgctggggacggcaccatcgtcttctgtgccgagcaggtaggagtggtgtgggtctacctccctgatgggagtcgctggagcaacccttc
ctggacctcaagaacatcgtgttgaccacccatggatcggggatgagagaggtcttgggtggtctttcaccccaaatccgccacaatcgcaagtctctattattatctgctggacaagaagaaggtagaaaagatccgaattagtgagatgaaggtttctcgggctgaccta
caagctgacctgaaatcagagagggtcatctggagattgaagaaccagcctcaaacataatggcggaactcttttggcctggtatgctatatgtacatatctactggggacggggacagctggagatcccttggcctgtttgaaatgctcagaacaaaagttccctgctggg
aaaagtfttaaggatcgtgtgaacagggcaggctcacatggcaagcgttacgagtcctcggacaatccattgttctgagccaggggccaccccgccatctatgcctatgggatcaggaaacatgtggcgttgtgtgtgtgaccgaggggaccccatcacgcgcagggccgag
gccgatattctgtgggacgtgggccaacaggttgaatg