

## ACTR3B-ARP3 actin-related protein 3 homolog B (yeast) Gene

Reference: PTXBC015207

|                |                                                                                                                                                                                                            |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ARP11; ARP3BETA; actin-related protein 3B; ARP3-beta; actin-like protein 3B; actin-related protein 3-beta; actin-related protein ARP4; actin-related protein Arp11; ARP3 actin related protein 3 homolog B |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ARP3 actin-related protein 3 homolog B (yeast)                                                    |
| <b>GENE ID</b>          | 57180                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015207                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ACTR3B                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcaggctccctgcctccctgcgtgggactgtggcaccgggtataccaagcttggctacgcaggcaacactgagccccagttcattatcttcattgtattgccatcagagagtcagcaaaagtagttgaccaagctcaaggagagtggtgaggggagttgatgaccttgacttttca  
tagggagatgaagccatcgataaacctacatafctacaaagtggccgatacgacatggaaatcattgaagactgggacttctatgaaagggtcatggagcaagtggtttttaaatactctcagctgaacctgaggaccattatttttaatagacagaacctccactcaatacaccagaaaacaga  
gagtatcttgcagaaattatgttgaatcatttaacgtaccaggactctacattgcagttcaggcagtgctggccttggcggcatcttgacatctcgacaagtggtggaacgtacgttaacggggatagtcattgacagcgagatggagtcacccatgttatccagtggcagaagggtatgt  
aattggaagctgcacatcaaacacatcccgattgcaggtagagatattacgtatttcattcaacagctgtaaggagagggagggtgggaatccctcctgagcagtcactggagaccgcaaaagccattaaggagaatactgttaccatttgccttgatagtcagggaatttccaagtatgat  
gtggatccccggagtgatcaaacagctacacgggtatcaatgcgatcaaccagaagaattgttatagacgttggttacgaaagattctgggacgtgaaatattcttaccggagtttgccaaccagactttatggagtcacatcagatgtgtgatgaagtaatacagaactgcccc  
atcgatgtgcggcgccgctgtataagccgagttcttcaggtctgccacaccaagaaggactatgaagagtagcgggccagcatctgccgccacaacccgctcttggagtcattgtcctag