



Reference: PTXBC015234

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARAT; GS1999FULL; HMFN1045; diacylglycerol O-acyltransferase 2; acyl-CoA retinol O-fatty-acyltransferase; diacylglycerol O-acyltransferase homolog 2; diacylglycerol O-acyltransferase-like protein 2; diglyceride acyltransferase 2		

For research use only.

GENE DESCRIPTION	diacylglycerol O-acyltransferase homolog 2 (mouse)
GENE ID	84649
ORIGIN SPECIES	Human
ACCESSIONS	BC015234
VECTOR	pDONR223
NCBI SYMBOL	DGAT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]