



Reference: PTXBC015244

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAD2B; POLZ2; REV7; mitotic spindle assembly checkpoint protein MAD2B; MAD2 (mitotic arrest deficient, yeast, homolog)-like 2; MAD2-like protein 2; REV7 homolog; hREV7; mitotic arrest deficient 2-like protein 2; mitotic arrest deficient homolog-like 2; po		

For research use only.

GENE DESCRIPTION	MAD2 mitotic arrest deficient-like 2 (yeast)
GENE ID	10459
ORIGIN SPECIES	Human
ACCESSIONS	BC015244
VECTOR	pDONR223
NCBI SYMBOL	MAD2L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]