

UBE2Q1-ubiquitin-conjugating enzyme E2Q family member 1 Gene

Reference: PTXBC015316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GTAP; NICE-5; PRO3094; UBE2Q; ubiquitin-conjugating enzyme E2 Q1; E2 ubiquitin-conjugating enzyme Q1; galactosyl transferase-associated protein; protein NICE-5; ubiquitin carrier protein Q1; ubiquitin conjugating enzyme E2Q family member 1; ubiquitin-conj		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2Q family member 1
GENE ID	55585
ORIGIN SPECIES	Human
ACCESSIONS	BC015316
VECTOR	pDONR223
NCBI SYMBOL	UBE2Q1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaactctcaccaacagggctggagcagtgctactccatagagtcagtgatcatgcagatcagtgccacactgtgaaggggaaagcacgagtcagtttgagccaacaatctcaatacagtcgtgacaagagcacagcagtcctacaagtccttggtgcagatccacgaaaaa
aacggctggtacacacccccaaaagaagacggctaa