

ST13-suppression of tumorigenicity 13 (colon carcinoma) (Hsp70 interacting protein) Gene

Reference: PTXBC015317

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AAG2; FAM10A1; FAM10A4; HIP; HOP; HSPABP; HSPABP1; P48; PRO0786; SNC6; hsc70-interacting protein; Hsp70-interacting protein; aging-associated protein 2; heat shock 70kD protein binding protein; progesterone receptor-associated p48 protein; renal carcinoma		

For research use only.

GENE DESCRIPTION	suppression of tumorigenicity 13 (colon carcinoma) (Hsp70 interacting protein)
GENE ID	6767
ORIGIN SPECIES	Human
ACCESSIONS	BC015317
VECTOR	pDONR223
NCBI SYMBOL	ST13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttttaaatcctttaaaacactcaccatataaaactgcatttgagcttggtgttcttttgtaatgtgtagagttctcttcgaaatgccagtggtgacttggttaactcaagaacagtttcttctggattcctatttgatttattaacctaattatattctaa