

PBRM1-polybromo 1 Gene

Reference: PTXBC015323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BAF180; PB1; protein polybromo-1; BRG1-associated factor 180; polybromo-1D; polybromo 1		

For research use only.

GENE DESCRIPTION	polybromo 1
GENE ID	55193
ORIGIN SPECIES	Human
ACCESSIONS	BC015323
VECTOR	pDONR223
NCBI SYMBOL	PBRM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagaagaagatagtgaggtcattgaacctcttctctacgtcagctcagacccccctggccagtgaagctcctcctacacacccccacagctacccccaaagtctgccaaggcagtgcaagaagggaaggtccaaacggaaatcaacatgagtggtctacatcgttcagcagtgagatgagggctgtgattaaggcccaacacccagactactcttctggggagctcagccgctgtgtggggacagaatggagaatcttgagacagccaagaagcagaatatgaaggcatgaggggtggtatccgcccaggtctccaccttgcagggccagttgatggcctgttagcatgggcagcatgcagccacttcacctgggggctccacccaccatcttccgcccaggtgtgcttggcctccgggcatccaccacgggtgtgatgaaccaaggagtgccctctatgttagggactccagcaccaggtggaagtcctatggacaacaggtgggagttttgggctccaggcagcagggcaccacctccatccgggcccacatccagctggacccccgtctacacagcagccaacaacacccatgttttagctccccacaaagaccagcggcctcttcaactcagaggcctacctaataacattgaaggactcagtcgggagtcacaagcattagcaagtgggatcagacactggcagctcgaagacgcgacgtccattgtcgaagaacaggagagccgctacacctctcactggctgaaaagcaaggggcccaaccaccatggcagatgccctctggcgcttcgagatttgatgtccgggacaccccaacattcgccaaacataacctaagaattgttaa