

SMR3B-submaxillary gland androgen regulated protein 3B Gene

Reference: PTXBC015327

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PBII; PRL3; PROL3; SMR1B; submaxillary gland androgen-regulated protein 3B; proline rich 3; proline-rich peptide P-B; proline-rich protein 3; salivary proline-rich protein; submaxillary gland androgen regulated protein 3 homolog B; submaxillary gland andr		

For research use only.

GENE DESCRIPTION	submaxillary gland androgen regulated protein 3B
GENE ID	10879
ORIGIN SPECIES	Human
ACCESSIONS	BC015327
VECTOR	pDONR223
NCBI SYMBOL	SMR3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatcactgacttggatcttgggcttgggctcttgcagcgtgtttcacacctggtagagagcaaaagagggcccgaccatccacctggaccgctggtctctctcaaccttttggccaggaattgttcaccacctctctccacccatggtccagggagaaatcccacctcc
tcctcccgccacctatggtccagggatatttccaccacccctctctcaacctaa