

TAGLN3-transgelin 3 Gene

Reference: PTXBC015329

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NP22; NP24; NP25; transgelin-3; neuronal protein 22; neuronal protein NP25; transgelin 3		

For research use only.

GENE DESCRIPTION	transgelin 3
GENE ID	29114
ORIGIN SPECIES	Human
ACCESSIONS	BC015329
VECTOR	pDONR223
NCBI SYMBOL	TAGLN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaacaggggcccagctatggcttaagccgagaggtgcaggagaagatcgagcagaagtatgatcgccacgtggagaacaagctgtggactggatcctgcagtgcgccgaggacatagagcaccgcccccgaggcccatcttcagaaatggtaatggacgg
gacggtcctgtgcaagctgataaatagttataccaccaggacaagagccatacccaagatctcagagtaaaagatggcttttaagcagatggagcaaatctccagttcctaaaagctgcggagacatggtgtcagaaccaccgacatcttcagacggtggatctatgggaaggga
aggacatggcagctgtgcagaggacctgatggctttaggcagcgtgtgcagtcaccaaggatgatggctctatcggggagagccatcctgtttcacaggaaagccagcagaatcgagaggcttttccgaggagcagcttgcaggggacagaacgtaataggcctgcagatggg
cagcaacaagggagcctccaggcggcatgacaggggtacgggtacggatgccagcagatcatgtag