

HOXB7-homeobox B7 Gene

Reference: PTXBC015345

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HHO.C1; HOX2; HOX2C; Hox-2.3; homeobox protein Hox-B7; homeo box 2C; homeo box B7; homeo box c1 protein; homeobox protein HHO.C1; homeobox protein Hox-2C; homeobox B7		

For research use only.

GENE DESCRIPTION	homeobox B7
GENE ID	3217
ORIGIN SPECIES	Human
ACCESSIONS	BC015345
VECTOR	pDONR223
NCBI SYMBOL	HOXB7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttcattgtattatgcgaatgctttattttctaaatccagcctcaagttcgttttctaccggagccttccagaacaactctgtgtgctttgttccaacccccagcgccgggtatggagcgggttcggcgcttcttcgccgctcgtatgcagggtgttaccggcgggg
ggggcatggcgccagagcgcgccggcgctctacgcggcggtatgggtcgcgagcggagttccttcaacatgcactgcgcgcctttgagcagaacctctccgggtgtgtccggcgactccgcaaggcgcgccggcgccaaggagcagagggactggactggcgccg
agagtaacttccgatctaccttggatgcgaagctcaggaactgaccgcaaacgagggccgacacctacacccgctaccagacctggagctggagaagaatttactacaatcgctacctgacgcggcgccggcgcatcgatcgccgacacgctctgctcacggaagac
agatcaagatttggttcagaaccggcgcatgaagtggaaaaaggagaacaagaccgcccggggaccaccggccaagacagggctgaagcagaggagggaaggaagagtga