

EBI3-Epstein-Barr virus induced 3 Gene

Reference: PTXBC015364

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-27B; IL27B; IL35B; interleukin-27 subunit beta; EBV-induced gene 3 protein; Epstein-Barr virus induced gene 3; IL-27 subunit beta; IL27 subunit; IL35 subunit; cytokine receptor; epstein-Barr virus-induced gene 3 protein; Epstein-Barr virus induced 3		

For research use only.

GENE DESCRIPTION	Epstein-Barr virus induced 3
GENE ID	10148
ORIGIN SPECIES	Human
ACCESSIONS	BC015364
VECTOR	pDONR223
NCBI SYMBOL	EBI3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccccgagcttctctgcccctgtcctctggccagctgcccgcctgagtggaaggaaaggccccagcagctctgacactgccccgggtgcaatgccgagcctctcggtaccgatgccgtggattgctcctggaccctgccgctgctccaaactccaccagccccgtgtcctctatggcaggtacaggtcgcctggtgctgcccggggccacagctggccctgctgcagcagacgccaacgtccaccagctgcacatcacggatgccagctgttctccatggctccctacgtgctcaatgtcaccgcccgtccacccctggggctccagcagcagcttcgtgcctttcataacagagcacatcatcaagcccacccctccagaaggcgtgctgctaaagccccctgctgagcgccagctacaggtgcagtgaggcctccgggtcctggccctccagagatcttctcactgaagtactggatccgttacaagcgtcaggagctgctgctccaccgggtggggccattgaagccagctctcctcctcagggctgtgctgccccgagccaggtactactgccaagtggcggctcaggacctcacagactacggggaactgagtgaactgagctctccccgccactgcccaatgagcctgggcaagtag