



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM99-transmembrane protein 99 Gene

Reference: PTXBC015365

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

transmembrane protein 99

For research use only.

GENE DESCRIPTION	transmembrane protein 99
GENE ID	147184
ORIGIN SPECIES	Human
ACCESSIONS	BC015365
VECTOR	pDONR223
NCBI SYMBOL	TMEM99
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]