

DPM2-dolichyl-phosphate mannosyltransferase polypeptide 2, regulatory subunit Gene

Reference: PTXBC015374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDG1U; dolichol phosphate-mannose biosynthesis regulatory protein; DPM synthase complex subunit; DPM synthase subunit 2; dolichol-phosphate mannose synthase subunit 2; dolichyl-phosphate mannosyltransferase polypeptide 2, regulatory subunit; dolichyl-phos		

For research use only.

GENE DESCRIPTION	dolichyl-phosphate mannosyltransferase polypeptide 2, regulatory subunit
GENE ID	8818
ORIGIN SPECIES	Human
ACCESSIONS	BC015374
VECTOR	pDONR223
NCBI SYMBOL	DPM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacggggacagaccaggtgtgggactcggcctcgtcgcgttagcctgatcatcttcactactacacgcctgggtgattctcttggtatgtcattctccccgtccgtgtcaccctccccgagccctggcaccgccagagcaactactataggtcttaggcacggcgtggc
ttcattgctgcctcatctcttga