

RPS5-ribosomal protein S5 Gene

Reference: PTXBC015405

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	40S ribosomal protein S5; ribosomal protein S5		
For research use only.			

GENE DESCRIPTION	ribosomal protein S5
GENE ID	6193
ORIGIN SPECIES	Human
ACCESSIONS	BC015405
VECTOR	pDONR223
NCBI SYMBOL	RPS5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgagtgaggagacagcaccagcggtggcagagacccagacatcaagctcttgggaagtggagcaccgatgatgtgcagatcaatgacattccctgcaggattacattgcagtgaaggagaagtatgccaagtacctgcctcacagtgcagggcggtatgccgcaaac
gctccgcaaagctcagtgctccattgtggagcgcctcaactaatgatgatgcacggccgcaacaacggcaagaagctcatgactgtgcgcacgtcaagcatgccttcgagatcatacaccctgctcacaggcgagaaccctctgcaggtcctgtgaacgccatcatcaacagtgg
tccccgggaggaactccacacgcattggcgccgggactgtgagacgacaggctgtggatgtgtccccctgcgccgtgtgaaccaggccatctggctgtgtgcacaggcgctcgtgaggtgccttcgggaacattaagaccattgctgagtgcttgccagatgagctcatcaatg
ctgccaagggtcctcgaactcctatgccattaagaagaaggacgagctggagcgtgtggccaagtccaaccctga