

IFI27L1-interferon, alpha-inducible protein 27-like 1 Gene

Reference: PTXBC015423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM14B; ISG12C; interferon alpha-inducible protein 27-like protein 1; ISG12(c) protein; family with sequence similarity 14, member B; interferon-stimulated gene 12c protein; interferon alpha inducible protein 27 like 1		

For research use only.

GENE DESCRIPTION	interferon, alpha-inducible protein 27-like 1
GENE ID	122509
ORIGIN SPECIES	Human
ACCESSIONS	BC015423
VECTOR	pDONR223
NCBI SYMBOL	IFI27L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaaggagagtgatgggactcaggcaggctgctgtacagctgtgtgctggaggagttgtgctgtgggactgtgctgtggcgctcagtgccatgggcttcacctcagtaggaatcgcgcacctcctcatagcagccaagatgatgtctacagcagccattgccaacggggcggagttgctgtggcagctgtgtgctattctgcagtcagtggggcagctggactctgtgacatctaaagttaacgggggctttgtgggacagctcttgggctgtgggttaccccccttcagctga