

MRPS33-mitochondrial ribosomal protein S33 Gene

Reference: PTXBC015462

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-139; MRP-S33; PTD003; S33mt; 28S ribosomal protein S33, mitochondrial; mitochondrial ribosomal protein S33		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S33
GENE ID	51650
ORIGIN SPECIES	Human
ACCESSIONS	BC015462
VECTOR	pDONR223
NCBI SYMBOL	MRPS33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctcccttcagaatatgccttcgcgtctcagtcgccggcgtatttggtgaagtcaccaggcctactaattccaagtctatgaaagtgtgaaactgtttagtgaactgcccttgccaagaagaaggagacttatgattggatccaaatcaccacacttacgtgaactcatgca
gacgtccgatttcttgactctacagagatgagcatcaggattttatgtagcgaacaaacgactaaagaagcttcgtggaagggagaaccaaagaaggagaagggaagagagcagcaaaagggaatatg