

C19orf53-chromosome 19 open reading frame 53 Gene

Reference: PTXBC015465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC023; LYDG10; leydig cell tumor 10 kDa protein homolog; chromosome 19 open reading frame 53		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 53
GENE ID	28974
ORIGIN SPECIES	Human
ACCESSIONS	BC015465
VECTOR	pDONR223
NCBI SYMBOL	C19orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaggggcagcgcaagtttcaggcgacaaacccgcaagagtaagacggcagcgagcctctgaaaagaatcgggcccaagaaaaggcggtcgtgttatcgtcccaagaaggcgcgctcgtgcagcagcaaaagctcaagaagaacctagaagtcggaatcgg
aagaagatcgaacatgacgtggtgatgaaagccagcagcagcctgcccagaagctggcactgctgaaggcccaagcaagaagaaggcgagctgccgccacctctccaagacaccttctga