



Reference: PTXBC015497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFTR-2; RTEF1; TCF13L1; TEF-3; TEF3; TEFR-1; hRTEF-1B; transcriptional enhancer factor TEF-3; TEA domain family member 4; related transcription enhancer factor 1B; transcription factor 13-like 1; transcription factor RTEF-1; transcriptional enhancer facto		

For research use only.

GENE DESCRIPTION	TEA domain family member 4
GENE ID	7004
ORIGIN SPECIES	Human
ACCESSIONS	BC015497
VECTOR	pDONR223
NCBI SYMBOL	TEAD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]