

## EDF1-endothelial differentiation-related factor 1 Gene

Reference: PTXBC015500

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFAP280; EDF-1; MBF1; endothelial differentiation-related factor 1; multiprotein bridging factor 1; endothelial differentiation related factor 1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	endothelial differentiation-related factor 1
<b>GENE ID</b>	8721
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC015500
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	EDF1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggccgagagcactgggacacgggtgacggctgctgcgcaagaaggccctacggccgcccagcccaaatccaagcaggctatcttagcggcacagagacgaggagaagatgtggagactccaagaatgggctgctggccagaacaacaacattctattaccaagaacacgg  
ccaagctggaccgggagacagaggagctgcacatgacagggtagcctggaggtgggcaaggatgacgcaaggtcggcagagcaagggcttacgagaaggacctggccacgaaaatcaatgagaagccacaggtgatcgcggactatgagagcggacggccataccc  
aataaccaggtgctggcaaaatcagcgggccattggcctcaagctcggggaaggacattggaagcccatcgagaagggcctagggcgaaatga