

AGR2-anterior gradient homolog 2 (Xenopus laevis) Gene

Reference: PTXBC015503

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AG2; GOB-4; HAG-2; HEL-S-116; PDIA17; XAG-2; anterior gradient protein 2 homolog; AG-2; HPC8; anterior gradient homolog 2; epididymis secretory protein Li 116; protein disulfide isomerase family A, member 17; secreted cement gland homolog; secreted cement		

For research use only.

GENE DESCRIPTION	anterior gradient homolog 2 (Xenopus laevis)
GENE ID	10551
ORIGIN SPECIES	Human
ACCESSIONS	BC015503
VECTOR	pDONR223
NCBI SYMBOL	AGR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaaattccagtgtagcattctgtgctctgtggccctctctacactctggccagagataccacagtcacaaacctggagcctctctgacccaaactgcccagaccctctccagaggttgggtgaccaaactcatctggactcagacatatgaagaagctctataataatccagacaagcaacaaccccttgatgattatcatcacttggaatgagtgcccacacagtcagctttaaagaagtggttctgctgaaataaagaatccagaattggcagagcagttgtctctcctcaatctggttatgaacaactgacaaacaccttctctgatggccagtatgtcccagagattatgtttgacccatctctgacagtttagagccgatatactggaagatatcaaacctctctatgcttacgaacctgcagatacagctctgttgcctgacaacatgaagaagctctcaagttgctgaagactgaattgtaa