

IL6-interleukin 6 (interferon, beta 2) Gene

Reference: PTXBC015511

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BSF-2; BSF2; CDF; HGF; HSF; IFN-beta-2; IFNB2; interleukin-6; B-cell differentiation factor; B-cell stimulatory factor 2; CTL differentiation factor; hybridoma growth factor; interferon beta-2; interleukin BSF-2; interleukin 6		

For research use only.

GENE DESCRIPTION	interleukin 6 (interferon, beta 2)
GENE ID	3569
ORIGIN SPECIES	Human
ACCESSIONS	BC015511
VECTOR	pDONR223
NCBI SYMBOL	IL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactccttccacaagcgcttcgggtccagttgccttccctggggctgctctggtgttgctgctgcttccctgccccaggtacccccaggagaagattccaaagatgtagccgccccacagacagccactcactcttcagaacgaattgacaaacaaatcggtacatcctg
acggcatctcagccctgagaaaggagacatgtaacaagagtaacatgtgtgaaagcagcaaaaggcactggcagaaaaaacctgaaccttccaaagatggctgaaaaagatggatgcttccaatctggattcaatgaggagacttgctggtgaaaaatcatcactggtcttttggagttt
gaggtatactagagtacctccagaacagatttgagagtagtgaggaaacagccagagctgtgcagatgagtagacaaaagtcctgatccagttcttcagaaaaaggcaaaagatctagatgcaataaccaccctgaccaaccacaaatgccagcctgctgacgaagctgcaggcaca
gaaccagtggtgcaggacatgacaactcatctcattctgagcagcttaaggagttctgcagtcagcctgagggctcttcggcaaatgtag