

PDCD5-programmed cell death 5 Gene

Reference: PTXBC015519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TFAR19; programmed cell death protein 5; TF-1 cell apoptosis-related protein 19; TF1 cell apoptosis-related gene 19; TFAR19 novel apoptosis-related		

For research use only.

GENE DESCRIPTION	programmed cell death 5
GENE ID	9141
ORIGIN SPECIES	Human
ACCESSIONS	BC015519
VECTOR	pDONR223
NCBI SYMBOL	PDCD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacgaggagcttgaggcgtgaggagacagagctggcggcgtgagcccaacacggggtcctggtgatggcccaacaggaagcaaacacagggaagcagaatgagaacagtatcttagcccaagttctggatcagtcggccggccagggttaagtaact
tagcactgttaagcctgaaaaactaaagcagtagagaattaccttatcacagatggcaagatattgacaactaagtgaagaaggtatcagaacaagggttaataagaatccttaaaaaagtaagccaacaacagaaaagacaacaacagtgaattcaacagaagaaagtaatgactct
gatgaagatgacgattattga