

LOC222699-transducer of ERBB2, 2 pseudogene Gene

Reference: PTXBC015543

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	transducer of ERBB2, 2 pseudogene
GENE ID	222699
ORIGIN SPECIES	Human
ACCESSIONS	BC015543
VECTOR	pDONR223
NCBI SYMBOL	LOC222699
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggaccccggtggggagctggccgctaagcgaagtggcctgacagtggagatgtgcgggccaatgtgcctgaggagctgagcatatggattgacccctcgggggtgctaccagattggtgagaaggagcagtgaaagtgtgtacctggatgacagtacggctgtgggg
ccccggagctggacatgaagatcaagagcagtttcaactctgacgaccagatgctctttctcctcggaagccaggacagctccttgccaactccccgtcgtatgctttgccatcaccagccccaccttcacctctgctogcccagcccaccttcacatggactcctttgctgcca
ccaaattag